

⚠ Instructions for Use – This leaflet contains important product use and safety information. Please read carefully and retain these instructions for future reference.

PRODUCT NAME:
ADVANTAGE® RADIOACTIVE SEED LOCALISATION
(I-125 RADIOACTIVE SEED LOCALISATION NEEDLE)

Description:
 The IsoAid ADVANTAGE RSL® I-125 Radioactive Seed Localisation Needle [RSLN] is a pre-sterilised 18-gauge stainless steel needle containing a low activity I-125 Iodine Seed (Advantage™ I-125 source). The ADVANTAGE™ I-125 source consists of a laser-welded Titanium capsule, containing Iodine-125 chemically affixed (adsorbed), as silver iodide, onto a silver rod which acts as an x-ray detectable marker. The needle tip is occluded with bone wax and the iodine seed is loaded loose or stranded and is provided with or without a trailing spacer. The stainless-steel needles are provided in 5cm, 7cm and 12cm lengths.

The ADVANTAGE RSL is intended for use by licensed and trained healthcare providers. It is available by prescription only.

Patients may be exposed to the following absorbable materials as part of the RSLN device: PGLA (90% glycolide/10% L-Lactide) polymer (if spacers are used). PLDLA 70L/30DL (if seeds are stranded) and bone wax which is comprised of beeswax, paraffin and isopropyl palmitate. All materials have been tested and comply with biocompatibility standards for implantable materials.

Indications for Use:
 The I-125 Radioactive Seed Localisation Needle is indicated to aid in the diagnosis of non-palpable tumours, lesions or associated lymph nodes in the breast by defining tumour, lesion or node location with a radioactive seed in preparation for excision. The Advantage RSL (RSLN) is intended to be used on adult individuals with non-palpable tumours/lesions/nodes. The radioactive seed is intended to be excised within thirty (30) days of implant.

Intended Use/Intended Purpose:
 The I-125 Radioactive Seed Localisation Needle is intended as a temporary implant to aid in localisation and excision of a tumour, lesion or associated lymph nodes in the breast. It is intended to be used with or without an absorbable strand and spacer.



Physical Characteristics:
 Iodine-125 has a half-life of 59.41 days and decays by electron capture with the emission of characteristic photons and electrons. The principal photon emissions are 27.2 KeV, 27.5

KeV, 31.0 KeV and 35.5 KeV with an average energy of 28.5 KeV. Table 1 shows the decay of I-125 seeds

Calibration:
 ADVANTAGE™ I-125 sources are calibrated by direct comparison against a standard source of the same model that has been calibrated by the US National Institute of Standards and Technology for Air Kerma Strength. The resulting calibration is reported in Air Kerma Strength (µGy m²/h) as well as Apparent Activity (mCi).

ADVANTAGE™ I-125 sources are calibrated to the NIST SK99std WAFAC standards for I-125 seeds.

Sterilisation/Single Use Only:
 The Radioactive Seed and Localisation Needle are sterilised with a Sterility Assurance Level of 10⁻⁶ by Ethylene Oxide gas. The sterile packaging has a one hundred and eighty (180) day shelf-life. If the product's expiry date has been exceeded, the product is considered non-sterile and therefore cannot be used. **Do not re-sterilise the product.**

Table 1. Iodine-125 Decay for RSLN

| Day | Decay Factor | Day | Decay Factor | Day | Decay Factor | Day | Decay Factor |
|-----|--------------|-----|--------------|-----|--------------|-----|--------------|
| 0 | 1.000 | 31 | 0.697 | 62 | 0.485 | 93 | 0.338 |
| 1 | 0.988 | 32 | 0.688 | 63 | 0.480 | 94 | 0.334 |
| 2 | 0.977 | 33 | 0.680 | 64 | 0.474 | 95 | 0.330 |
| 3 | 0.966 | 34 | 0.673 | 65 | 0.469 | 96 | 0.326 |
| 4 | 0.954 | 35 | 0.665 | 66 | 0.463 | 97 | 0.323 |
| 5 | 0.943 | 36 | 0.657 | 67 | 0.458 | 98 | 0.319 |
| 6 | 0.932 | 37 | 0.649 | 68 | 0.452 | 99 | 0.315 |
| 7 | 0.922 | 38 | 0.642 | 69 | 0.447 | 100 | 0.311 |
| 8 | 0.911 | 39 | 0.634 | 70 | 0.442 | 101 | 0.308 |
| 9 | 0.900 | 40 | 0.627 | 71 | 0.437 | 102 | 0.304 |
| 10 | 0.890 | 41 | 0.620 | 72 | 0.432 | 103 | 0.301 |
| 11 | 0.880 | 42 | 0.613 | 73 | 0.427 | 104 | 0.297 |
| 12 | 0.869 | 43 | 0.606 | 74 | 0.422 | 105 | 0.294 |
| 13 | 0.859 | 44 | 0.599 | 75 | 0.417 | 106 | 0.290 |
| 14 | 0.849 | 45 | 0.592 | 76 | 0.412 | 107 | 0.287 |
| 15 | 0.839 | 46 | 0.585 | 77 | 0.407 | 108 | 0.284 |
| 16 | 0.830 | 47 | 0.578 | 78 | 0.403 | 109 | 0.280 |
| 17 | 0.820 | 48 | 0.571 | 79 | 0.398 | 110 | 0.277 |
| 18 | 0.811 | 49 | 0.565 | 80 | 0.393 | 111 | 0.274 |
| 19 | 0.801 | 50 | 0.558 | 81 | 0.389 | 112 | 0.271 |
| 20 | 0.792 | 51 | 0.552 | 82 | 0.384 | 113 | 0.268 |
| 21 | 0.783 | 52 | 0.545 | 83 | 0.380 | 114 | 0.265 |
| 22 | 0.774 | 53 | 0.539 | 84 | 0.375 | 115 | 0.261 |
| 23 | 0.765 | 54 | 0.533 | 85 | 0.371 | 116 | 0.258 |
| 24 | 0.756 | 55 | 0.526 | 86 | 0.367 | 117 | 0.255 |
| 25 | 0.747 | 56 | 0.520 | 87 | 0.362 | 118 | 0.252 |
| 26 | 0.738 | 57 | 0.514 | 88 | 0.358 | 119 | 0.250 |
| 27 | 0.730 | 58 | 0.508 | 89 | 0.354 | 120 | 0.247 |
| 28 | 0.721 | 59 | 0.502 | 90 | 0.350 | | |
| 29 | 0.713 | 60 | 0.497 | 91 | 0.346 | | |
| 30 | 0.705 | 61 | 0.491 | 92 | 0.342 | | |

In Vivo Characteristics:
 During the excision procedure, the seed provides a radioactive localisation point and acts as a marker to aid in the location and excision of the lesion. Verify the removal of seed at time of excision of the tumour/lesion/node using a gamma probe or similar instrument intended for isotope detection.

Instructions for Safe Use:
 The radioactive seed is introduced via an 18-gauge needle using standard ultrasound or radiography guidance. Once guided to the desired location of the lesion, the seed is deployed through the bone wax with the aid of the needle stylet.

If multiple lesions utilise more than one seed, then each seed shall be a minimum of 2 cm apart. Ultrasound or radiography confirms the appropriate placement of the seed.

The seed is intended to be removed during the excision procedure.

For single use only. Do not re-use. Re-using this medical device bears the risk of cross-patient contamination, as medical devices – particularly those with long and small lumina, joints or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with micro-organisms which may lead to infectious complications.

Do not resterilise. After resterilisation, the sterility of the product is not guaranteed because of the indeterminable degree of microbiological contamination which may lead to infectious complications. Cleaning, reprocessing and/or resterilisation of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.

Radiation Protection & Handling:
 The 27- 35.5 KeV photons of I-125 are substantially absorbed by any high Z material but exhibit desirable penetration in tissue.

Half Value Layer Lead = 0.025 mm
 Half Value Layer Tissue = 20.0 mm

Exposure can be reduced by 99.9% with a thin sheet of lead (0.25 mm or 0.01 inch). The shielding of I-125 results in a reduction of exposure to attending medical personnel and visitors. I-125 sources should be handled only by those individuals trained by an authorising governmental agency in the safe use and handling of radio-isotopes.

- Direct contact with I-125 sources should be avoided. The use of vacuum or reverse action tweezers is recommended. Proper precautions must be taken when handling the sources.
- Personnel monitoring is required. Dosimetry monitors, such as TLD devices, should be used to monitor hand and whole body exposure. During preparation and source implantation procedures, all practical steps should be taken to keep exposure as low as reasonably achievable. Limited exposure time, increasing distance, careful planning of the administration procedure and use of shielded barriers should be considered in meeting this goal.

Accidental Damage:

Do not use the product if there is suspicion that the product is damaged or if the sterile barrier has been breached. It is possible through rough handling (abrasion, incision, etc.), high temperatures or crushing that a seed could rupture and leak. The internal components of the seed are non-toxic, but the area should be closed off immediately and personnel limited to avoid radioactive contamination. The damaged seeds should be placed in a sealed container and the area should be decontaminated. In accordance with radiation regulations, only authorised, specialised staff trained in handling radioactive substances may handle the I-125 seeds.

Accountability & Disposal:

Records of receipt, storage and disposal of Advantage™ I-125 sources should be maintained in accordance with government regulatory policies. I-125 sources should be strictly controlled and stored in a secured area.

When disposal is indicated, the Advantage™ I-125 sources should be transferred to an authorised radioactive waste disposal agency or returned to IsoAid for disposal. NOTE: IsoAid does not accept I-125 sources for return from the European Union. Advantage™ I-125 sources should not be disposed of in normal waste. Any discrepancies must be reported immediately to IsoAid Customer Service.

Licensing:

USA – State/Federal:

⚠ CAUTION: Federal (USA) and State law(s) restrict this device to sale by or on the order of a physician.

The Florida Department of Health (FDOH), Bureau of Radiation Control has approved this sealed source for distribution to persons licensed pursuant to Florida Administrative Code Chapter 64E-5, "Control of Radiation Hazard Regulations," Part VI or under equivalent licenses of the USNRC or issued by an Agreement State. IsoAid requires proof of USNRC radioactive materials license or respective government license as well as agreement state and licensing state information. Orders cannot be processed without license verification. Compliance with the applicable local, state, country, and/or government regulations concerning procurement, possession, use and disposal of radioactive materials is the responsibility of the customer.

Canada - Canadian Nuclear Safety Commission

⚠ CAUTION: Canadian National and Regional/Province law(s) restrict this device to sale by or on the order of a physician.

Application of REGDOC-2.12.3, Security of Nuclear Substances: Sealed Sources for typical uses of sealed

sources, Brachytherapy - low dose rate is a Category 4 source. Category 4 sources are very unlikely to permanently injure anyone. However, this amount of unshielded radioactive material, if not safely managed or securely protected, could possibly – although it is unlikely – temporarily injure someone who handled it or was otherwise in contact with it, or who was close to it for a period of many weeks. This Code of Conduct on the Safety and Security of Radioactive Sources was approved by the Board of Governors of the International Atomic Energy Agency (IAEA) on 8 September 2003. It replaces the version published (with the symbol IAEA/CODEOC/2001) by the IAEA in March 2001. It reflects the important findings produced by the International Conference on Security of Radioactive Sources held in Vienna in March 2003 (the Hofburg Conference). Member States to be encouraged to join and effectively implement these Conventions. Canada is already a signatory to these conventions, together with codes of conduct on non-proliferation, research reactors and the safety and security of radioactive sealed sources, along with the Comprehensive Nuclear Test-Ban Treaty.

Canadian Nuclear Safety Commission
280 Slater Street P.O. Box 1046
Station B Ottawa, Ontario K1P 5S9 CANADA
Tel.: 613-995-5894 or 1-800-668-5284 (in Canada only)
Facsimile: 613-995-5086 Email: info@cncs-ccsn.gc.ca
Web site: nuclearsafety.gc.ca

Australia - Australian Radiation Protection and Nuclear Safety Agency

⚠ CAUTION: Australian National and Regional/State law(s) restrict this device to sale by or on the order of a physician.

The establishment of a NRWMF is governed by the National Radioactive Waste Management Act 2012. A NRWMF also needs to adhere to the Environment Protection and Biodiversity Conservation Act 1999, the Nuclear Non-Proliferation (Safeguards) Act 1987 and the Australian Radiation Protection and Nuclear Safety Act 1998.

The proposed National Radioactive Waste Management Facility would be a controlled facility under the Australian Radiation Protection and Nuclear Safety Act 1998 (the Act).

Under the Act, licences are required to prepare a site for, construct, or operate a controlled facility. The decision to submit a licence application is a matter for the applicant.

Before an application is made to the CEO of ARPANSA for a licence to prepare a site for the National Radioactive Waste Management Facility, the applicant will have to obtain approval from the Minister for the Environment under the Environment Protection and Biodiversity Conservation Act 1999. Before any radioactive material is allowed to be transported it must be packed, shielded, labelled and marked as set out in the ARPANSA Code: Safe Transport of Radioactive Materials. This code is based on the International Atomic Energy

Agency's (IAEA) Regulations for Safe Transport of Radioactive Material. [nrwfmfsupport@arpansa.gov.au](mailto:nrwmfsupport@arpansa.gov.au); www.arpansa.gov.au

A radioisotope is considered to be for medical use when it is intended to be:

1. administered to humans or used for any therapeutic procedure or purpose in any planned exposure of humans to ionising radiation
2. used in any in vitro medical diagnosis or test
3. used in research which is either directly or indirectly related towards medical diagnosis or therapy in humans.

Note: Sealed and unsealed radioactive sources that are used to calibrate instruments in medical practices and pathology laboratories are also included as medical radioisotopes for permit purposes. The applicant/"end user" declares that he/she holds an appropriate licence issued by the relevant Commonwealth, State or Territory radiation regulatory authority to deal with the above radioisotopes. The applicant/"end user" also undertakes not to supply any of the above radioisotopes to an unapproved user. The applicant/"end user" should contact the relevant Commonwealth, State or Territory radiation regulatory authority for advice on legislative requirements.

medicalpermits@arpansa.gov.au; www.arpansa.gov.au.

ARPANSA, like other regulatory bodies in Australia and abroad, has been working on developing capability in holistic safety. Charged with the function of protecting the health and safety of people under the Australian Radiation Protection and Nuclear Safety Act 1998 (the Act), ARPANSA proposes to use a holistic approach to assess and monitor the safety of licence holders and applicants. These guidelines outline ARPANSA's vision and expectations for holistic safety.

Leak Testing:

ADVANTAGE I-125 Brachytherapy sources are 100% leak tested prior to shipment and have passed a leak test showing less than 185 Bq (5 nCi) of removable I-125 surface contamination as required by ISO 9978 "Radiation protection – Sealed radioactive sources". Advantage I-125 seeds do not require any additional leak testing provided the seeds are used within the use-by date.

Localisation Dosage and Administration:

The most commonly used source activity levels for localisation is between 0.1 mCi and 0.3 mCi.

⚠ Adverse Reactions:

- Any adverse reaction associated with tissue radiation damage may be associated with use of I-125 sources. Proper precautions must be taken when handling the sources.

- As with any surgical procedure, complications may occur, including: bruising, discomfort, prolonged bleeding or infection near the implant site.
- Although the risk of source migration is minimal, it can be significantly reduced through the use of stranding that links the seed and spacer together prior to implantation.
- Iodine Allergy
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the jurisdiction in which the user and/or patient is established.

⚠ Precautions:

- ⚠ Product should remain in leaded pouch until ready for use. Handle lead pouch and contents with care to prevent damage to product.
- ⚠ Use caution when patients are diagnosed with non-cancerous tumours/lesions. Implant and removal should occur within 24 hours to limit radiation exposure.

⚠ Contraindications:

- ⚠ Do not use radioactive seed localisation needles in neurological or cardiovascular tissues.
- ⚠ The RSLN is sold sterile. Use of a non-sterile device may compromise patient care. Do not re-sterilise.
- ⚠ Do not use a damaged seed or a seed that may have become damaged when using the device.
- ⚠ Do not come in direct contact with the I-125 source. Use vacuum or reverse action tweezers to handle the I-125 sources.
- ⚠ The needle is not to be used in an MRI environment.

⚠ Warnings:

- ⚠ Dispose of radioactive material per nuclear regulatory guidelines (for USA, 10 CFR 35.1000; for EU per EURATOM 1493/93)
- ⚠ Loss of a radioactive seed must be avoided. Protocols must be in place to ensure tracking of the seed throughout the process.
- ⚠ Any attempt to cut or segment stranded product may adversely result in radioactive contamination. Use product as intended.
- ⚠ Do not use if damaged. Discard if damaged during use or after use in accordance with waste disposal procedures.
- ⚠ Do not use when patients are pregnant or breast-feeding. An alternative non-radioactive device should be used to avoid radiation exposure.
- ⚠ Do not use on patients less than 18 years old. This product is intended for use in adults.
- ⚠ Do not use if needle is bent or broken.

- ⚠ Excessive Force is not required to expel seed.
- ⚠ Do not store without adequate leaded shielding /packaging
- ⚠ Healthy tissue may be exposed to the RSLN device during implantation and excision.

⚠ MR Conditional

The I-125 seed has been evaluated for safety in the MRI environment. The seeds are MR-Conditional as defined in ASTM F2503. The seeds have been tested for heating, migration and image artifact in the MRI environment. IsoAid seeds are made with a titanium shell with non-magnetic internal materials. Patients with the seeds may safely undergo MRI under the following conditions: 1) Static field of 3 T or less 2) Whole body SAR of 4 W/kg or less and head SAR of 3.2 W/kg or less 3) Normal or first level controlled mode of the MRI system for both RF and gradients 4) Maximum spatial gradient in the static field of 30 T/m (3000 Gauss/cm) 5) Maximum slew rate of the time-varying magnetic gradient for the seed is 200 [T/m/s], which is the high-end gradient slew rate and is worst-case for the seed that does not have any magnetic or transistors in the seed components, no conceivable negative impact.


⚠ The presence of other implants or the health state of the patient may require reduction of the MR limits.

⚠ Temperature rise of tissues surrounding the seed was calculated under a worst-case situation to be less than 50% above the background rise with no implant. Magnetic force and torque during MRI will be less than the values exerted by gravity. Image artifact is expected to extend less than 5 mm beyond the seeds.

⚠ CAUTION: Federal (USA) and State law(s) restrict this device to sale by or on the order of a physician.
















⚠ Use and Distribution in the EU is governed by EURATOM 2013/59 and 1493/93.

The RSLN Summary of Safety and Clinical performance is available on the European database on medical devices (EUDAMED), where it is linked to the Basic UDI: M936IRLS. <https://ec.europa.eu/tools/eudamed>.

 **The excised RSLN seed is considered biohazardous and must be contained and disposed of in accordance with universal precautions.**

The RSLN Product may be configured with or without a Spacer and/or Strand; and may be supplied in a 5cm, 7cm or 12cm stainless steel needle [where X = length of needle].

| Configurations (where "X" = cm in length, e.g. 5 cm, 7 cm, 12 cm) | Product Code |
|--|--------------|
| Stranded, no spacer | RSLN-X-SS |
| Stranded with spacer | RSLN-X-SS/S |
| Loose Load, no spacer | RSLN-X-LL |
| Loose Load with spacer | RSLN-X-LL/S |

| LEGEND | |
|---|------------------------------------|
| Caution: Consult Accompanying Documents | |
|  | Do Not Reuse |
|  | Consult Instructions for Use |
|  | Ethylene Oxide Sterilisation |
|  | Use by Date |
|  | Catalogue Number |
|  | Do not Resterilise |
|  | Biohazard |
|  | Radioactive |
|  | MR Conditional |
|  | Do not use if package is damaged |
|  | Importer |
|  | Medical Device |
|  | Manufacturer |
|  | Keep Dry |
|  | European Authorized Representative |

Clinical Benefit of RSL device for breast cancer patients:

1. Increased surgical accuracy: RSL allows surgeons to precisely locate and remove small, non-palpable breast tumours. This precision helps ensure that all cancerous tissue is removed while preserving as much healthy tissue as possible.
2. Patient comfort: Unlike traditional wire localisation, which involves a wire sticking out of the breast for several hours before surgery, RSL uses a tiny radioactive seed that is less uncomfortable for patients. The seed is placed using a needle under local anaesthesia and patients can go about their day without the discomfort of a protruding wire.
3. Reduced need for re-operation: By improving the accuracy of tumour removal, RSL reduces the likelihood of needing a second surgery to remove any remaining cancerous tissue.
4. Minimally invasive: The procedure is minimally invasive and can be performed under local anaesthesia, making it a safer option with fewer complications.

CE 2797

EU REP

Authorised Rep:

AJW Technology Consulting GmbH
Breite Straße 3
40213 Düsseldorf (Germany)
Tel: +49 211 54059 6030

Manufacturer:

IsoAid LLC

7824 Clark Moody Blvd
Port Richey, Florida 34668
United States of America
Ph: +1-727-815-3262